# Guidelines for cleaning, disinfection, and sterilization of the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor

This manual describes the procedures for cleaning, disinfecting, and sterilizing the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor used with the ultrasound transducers for Canon Medical Systems diagnostic ultrasound systems. For the precautions and detailed information concerning use of the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor, refer to the corresponding

operation manuals.

#### Trademarks

This manual may include trademarks or registered trademarks of companies other than Canon Medical Systems.

# 1. Safety Precautions

#### **1.1 Meaning of Signal Words**

In this manual, the signal words **DANGER**, **WARNING**, **CAUTION**, and *NOTICE* are used regarding safety and other important instructions. The signal words and their meanings are defined as follows. Please understand their meanings clearly before reading this manual.

Signal word	Meaning		
<b>ADANGER</b>	Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.		
<b>AWARNING</b>	Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.		
<b>▲CAUTION</b>	Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.		
NOTICE	Indicates a potentially hazardous situation which, if not avoided, may result in property damage.		

## **1.2 Safety Precautions**

Observe the following precautions to ensure the safety of patients as well as operators when performing cleaning, disinfection, or sterilization of the biopsy adaptor, water bag, coupler, or magnetic sensor securing adaptor.

<b>AWARNING</b>	Clean and sterilize the biopsy adaptor before and after each biopsy procedure. Note that local regulations may require that the biopsy adaptor be disinfected before sterilization. If the biopsy adaptor has not been properly sterilized, the biopsy adaptor may become a source of infection.		
	1. Precautions concerning cleaning, disinfection, and sterilization.		
	1) Observe the following precautions to prevent infection.		
	<ul> <li>Wear protective gloves when performing cleaning.*</li> </ul>		
	<ul> <li>Wear sterile protective gloves when performing disinfection or sterilization.*</li> </ul>		
	<ul> <li>Wear new protective gloves each time cleaning, disinfection, and sterilization are performed.</li> </ul>		
	<ul> <li>Clean the water bag and coupler and then disinfect or sterilize them before and after each examination.</li> </ul>		
	* Refer to the FDA's March 29, 1991 Medical Alert on Latex Products.		

AUTION	<ul> <li>2) After chemical cleaning, rinse the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor thoroughly with purified water. After chemical disinfection, rinse these items thoroughly with sterile water or deionized water.</li> <li>Chemical residues may be harmful to the human body.</li> </ul>
	3) After chemical cleaning or chemical disinfection, thoroughly dry the surfaces of the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor.
	4) After gas sterilization, perform aeration to degas the biopsy adaptor, water bag, and coupler. Gas residues on these items may be harmful to the human body.
	5) The efficacy of the cleaning solutions, disinfectants, and sterilizing gases is not guaranteed by Canon Medical Systems. Contact the manufacturers for information on the activity of the products.

<b>▲CAUTION</b>	6) To ensure the prevention of infection, confirm the effectiveness of each chemical for cleaning, disinfection, or sterilization based on the criteria (such as effective period, number of times of use, discoloration, and results of using the effectiveness test kit) described in the documentation provided by the relevant manufacturer.
	7) Do not clean, disinfect, or sterilize the biopsy adaptor, water bag, or coupler using chemicals or methods other than those specified in this guide. If chemicals or methods other than those specified in this guide are used, the item may not be properly cleaned, disinfected, or sterilized or may be damaged.
	8) The cleaning, disinfection, or sterilization conditions, such as the temperature and pressure, differ depending on the product. In addition, some products cannot be subjected to disinfection or sterilization procedures. Confirm the detailed conditions by referring to the tables in section 3.

## **1.3 Chemical Hazard (for USA only)**

Observe the following instruction in order to protect patients and operators from inflammation or poisoning by chemical substances.

<b>AWARNING</b> :	This product can expose you to chemicals including lead and lead compounds, which are known to the State of California to cause cancer, and phthalates, which are known to the State of California to cause birth defects or other reproductive harm. For more information go to
	<u>www.p65warnings.ca.gov</u> .

### **1.4 Precautions Concerning Restriction of Sale or Use** (for USA only)

United States law restricts this device to sale or use by, or on the order of a physician.

# 2. Cleaning, Disinfection, Sterilization

This section describes the methods for cleaning, disinfection, and sterilization.

## 2.1 Cleaning

- < <Items to be used: Protective gloves, cleaning solution, purified water, clean soft cloth or gauze, single-use sponge\*, soft singleuse brush>>
- \* The single-use sponge must not include any abrasive parts or contain any abrasive cleanser.
- Wear protective gloves to prevent infection.
   Wear new protective gloves each time cleaning is performed.

- (2) Confirm that all parts are disassembled. Disassemble all parts that are not yet disassembled. For further information on parts that can be disassembled, refer to the operation manuals for the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor.
- (3) Wash off all organic materials (such as blood or other bodily fluids) from the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor with purified water. A single-use sponge can be used to wash the water bag. The water bag is made of a very thin membrane. Do not apply excessive force to the water bag during washing. For the biopsy adaptor and coupler, a single-use sponge or soft single-use brush can be used. Do not use a hard brush because it may damage the biopsy adaptor or coupler. Use a soft single-use brush to wash the inside of the needle guide of the biopsy adaptor.

Use a new sponge or brush each time washing is performed. If a sponge or brush is reused, it may become a source of infection.

- (4) In accordance with the tables in section 3, immerse the biopsy adaptor, water bag, coupler, or magnetic sensor securing adaptor in a cleaning solution to dissolve or remove all remaining organic materials. Use a single-use sponge or soft single-use brush if necessary. The cleaning solution must reach the inside of the needle guide. If dried organic materials are present on an item, immerse it in the cleaning solution for a prolonged period.
- (5) Remove all residual organic materials and chemical from the biopsy adaptor, water bag, coupler, or magnetic sensor securing adaptor by rinsing with purified water.

- (6) Dry the surface of the biopsy adaptor, water bag, coupler, or magnetic sensor securing adaptor with clean soft cloth or gauze or an air blower. Be sure to dry the inside of the needle guide. Do not use heat for drying. Confirm that all organic materials and chemical have been completely removed.
- (7) For the detailed handling procedures for a chemical, refer to the documentation provided by the manufacturer of the chemical. In order to maintain the effectiveness of the chemical, be sure to observe the specified concentration, replacement interval, working temperature, and other requirements. Determine the effectiveness of the chemical based on the criteria (effective period, number of times of use, discoloration of the chemical, results of testing using effectiveness test kits, etc.) described in the documentation provided by the manufacturer of the chemical.
- (8) Confirm that the biopsy adaptor, coupler, and magnetic sensor securing adaptor are free of damage, deformation, or peeling and that they are functioning normally.
   Confirm that the water bag is free of damage, deformation, or peeling.

#### 2.2 Disinfection

Before disinfection, the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor must be cleaned. Note that local regulations may require that the biopsy adaptor, water bag, and coupler be disinfected before sterilization.

<<Items to be used: Sterile protective gloves, disinfectant, sterile water or deionized water, sterile soft cloth or gauze, soft single-use brush>>

- (1) Wear sterile protective gloves to prevent infection.Wear new sterile protective gloves each time disinfection is performed.
- (2) Disinfect the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor using the chemicals listed in section 3 "List of Chemicals". The disinfectants must reach the inside of the needle guide. Use a soft single-use brush if necessary.
  - (a) After performing disinfection, remove all residual chemical from the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor with sterile or deionized water. Do not reuse the sterile or deionized water.
  - (b) Dry the surfaces of the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor with sterile soft cloth or gauze or an air blower. Be sure to dry the inside of the needle guide. Do not use heat for drying.

Confirm that all organic materials and chemical have been completely removed.

For the detailed handling procedures for a chemical, refer to the documentation provided by the manufacturer of the chemical.

(3) To maintain the effectiveness of the disinfectant, ensure that the concentration, temperature, and other conditions specified in the documentation provided by the manufacturer are met.

To confirm the effectiveness of the disinfectant, use the criteria (such as effective period, number of times of use, discoloration, and results of using the effectiveness test kit) described in the documentation provided by the manufacturer.

(4) Confirm that the biopsy adaptor, coupler, and magnetic sensor securing adaptor are free of damage, deformation, or peeling and that they are functioning normally.
 Confirm that the water bag is free of damage, deformation, or peeling.

#### 2.3 Sterilization

Before sterilization, the biopsy adaptor, water bag, and coupler must be cleaned.

Note that local regulations may require that the biopsy adaptor, water bag, and coupler be disinfected before sterilization.

<< Items to be used: Sterile protective gloves, sterilant>>

- \* Some types of biopsy adaptor, water bag, or coupler cannot be sterilized, or the sterilization conditions may differ.
- Wear sterile protective gloves to prevent infection.
   Wear new sterile protective gloves each time sterilization is performed.
- (2) Drain the water from the water bag before sterilizing the water bag.
- (3) Sterilize the biopsy adaptor, water bag, or coupler using the chemicals listed in section 3 "List of Chemicals".
  - \* Place the biopsy adaptor, water bag, or coupler in a sterilization packing case, and then place it in the sterilizer.
- (4) In the case of high-pressure steam sterilization, heat the steam by performing prevacuuming three times and then start sterilization (sterilization temperature: 134°C (273.2°F), holding time: 18 minutes). Dry the items after sterilization.

- (5) After gas sterilization, perform aeration to remove all gas residues on the sterilized items. (Aeration is not required after hydrogen peroxide plasma sterilization.)
- (6) Confirm that the biopsy adaptor and coupler are free of damage, deformation, or peeling and that they are functioning normally.
   Confirm that the water bag is free of damage, deformation, or peeling.

# 3. List of Chemicals

This section lists the chemicals used when cleaning, disinfecting, or sterilizing the biopsy adaptor, water bag, coupler, and magnetic sensor securing adaptor. Refer to the lists for the available chemicals. In addition, refer to the manuals for the chemicals for the handling details.

NOTE For the transducers that support each item, refer to the operation manual supplied with the corresponding item.

# Cleaning

Chemical name/type Enzyme					
Country of origin :			3M™ Rapid		
Any: Any country		Trade name	CIDEZYME <sup>®</sup>	Multi-Enzyme	
FRA: France	e			cleaner	
DEU: Germ	any				
	d Kingdom	Manufacturer	1&1	3M Company	
USA: Unite	-				
JPN: Japan		Country of origin	Any	Any	
AUS: Austr		Concentration	0.8%	1%	
7105.71051	ana	(Dilution ratio)	(125)	(100)	
		Time	1 min.	10 min.	
		Temperature		mperature	
		Humidity		humidity	
		Pressure	Norma	pressure	
Water bag/	M	odel name <sup>*1</sup>			
coupler	UAWB-021A		ОК		
coupier	UAWB-022A				
	UAGL-001AHA				
	UAGL-002AHA	UAGV-027A			
	UAGL004A	UAGV-028A			
	UAGL-004AHA		OK		
Biopov	UAGV020A	UAGV-030A UAGV-034A			
Biopsy adaptor <sup>*2</sup>	UAGV-023A	UAGV-034A UAGV-037A			
adaptor	UAGV-035A				
	UAGV-036A			ЭК	
		UAGV-031A	ОК		
	UAGL-015A	UAGV-032A			
		UAGV-033A			
	UAFS-001A	UAFS-002A			
Magnetic	UAFS-003A	UAFS-004A	ОК		
sensor	UAFS-005A	UAFS-006A			
securing	UAFS-007A	UAFS-008A			
adaptor	UAFS-009A UAFS-011A	UAFS-010A UIFR-A500			
	UAL2-0114	UIFK-ADUU			

OK : Use of the chemical is permitted. \_\_\_\_ : Use of the chemical is not permitted.

Enzyme						
ANIOSYME DD1 <sup>*3</sup>	neodisher <sup>®</sup> mediclean <sup>*3</sup>	Tristel Pre-Clean Wipes <sup>*3</sup>	MetriZyme <sup>®</sup>	Klenzyme <sup>®</sup>	Instru- Zyme <sup>®</sup> Gel	
Laboratories ANIOS	DR. WEIGERT GmbH & Co. KG	Tristel Solutions Limited	Metrex Research, Inc.	STERIS Corporation	Summit Medical	
FRA	DEU	GBR	USA	USA	USA	
0.5% (200)	1% (100)	Working solutions	0.77% (130)	0.8% (125)	Working solution	
5 min.	10 min.	Wiping only	5 min.	5 min.	Wiping only	
	mperature		nperature		nperature	
	humidity		humidity		humidity	
Normal	pressure	Normal	pressure	Normal	Normal pressure	
ОК		ОК	ОК	ОК	ОК	
ОК		_	ОК	OK	_	
C	ОК		ОК	ОК	_	
ОК		_	ОК	_	_	
ОК		_	_	_	_	

- \*1: Note that some of the listed models may not be available in some regions. For the transducers that support each item, refer to the operation manual supplied with the corresponding item.
- \*2: Sterilize the biopsy adaptor before and after use.
- \*3: Not for use in the USA.

### Low/Middle-Level Disinfection -1 (Except for USA)

\* Note that local regulations may require that the biopsy adaptor, water bag, and coupler be disinfected before sterilization.

Country of Any: Any c FRA: Franc	ountry	Chemical name/type	Didecyl dimethyl ammonium chlorides/ Quaternary ammonium compounds	
DEU: Germ	-	Trade name	CLEANISEPT <sup>®</sup> WIPES	
USA: Unite JPN: Japan	d States	Manufacturer	Dr. Schumacher GmbH	
AUS: Austr		Country of origin	DEU	
AUS. AUST	alla	Concentration (Dilution ratio)	Working solution	
		Time	Wiping only	
		Temperature	Room temperature	
		Humidity	Normal humidity	
		Pressure	Normal pressure	
Water bag/	Mod	del name <sup>*1</sup>	ОК	
coupler	UAWB-021A UAWB-022A			
	UAGL-001AHA UAGL-002AHA UAGL004A UAGL-004AHA	UAGV-024A UAGV-027A UAGV-028A UAGV-029A	_	
Biopsy adaptor <sup>*2</sup>	UAGV020A UAGV-023A	UAGV-030A UAGV-034A UAGV-037A		
	UAGV-035A UAGV-036A		_	
	UAGL-015A	UAGV-031A UAGV-032A UAGV-033A	_	
Magnetic sensor securing adaptor	UAFS-001A UAFS-003A UAFS-005A UAFS-007A UAFS-009A UAFS-011A	UAFS-002A UAFS-004A UAFS-006A UAFS-008A UAFS-010A UIFR-A500	_	

OK : Use of the chemical is permitted. \_\_\_\_ : Use of the chemical is not permitted.

Isopropyl alcohol/ Didecyl dimethyl ammonium chlorides	Glyoxal/Glutaral/ Didecyl dimethyl ammonium chlorides		
WIP'ANIOS	Lysoformin <sup>®</sup> 3000		
Laboratories ANIOS	Lysoform Disinfection Ltd.		
FRA	DEU		
Working solution	2% (50)		
Wiping only	15 min.		
Room temperature	Room temperature		
Normal humidity	Normal humidity		
Normal pressure	Normal pressure		
ОК	ОК		
_	_		
_	_		
_	_		
_	_		

- \*1: Note that some of the listed models may not be available in some regions. For the transducers that support each item, refer to the operation manual supplied with the corresponding item.
- \*2: Sterilize the biopsy adaptor before and after use.

### Low/Middle-Level Disinfection -2 (High-level disinfection is also required in the USA.)

\* Note that local regulations may require that the biopsy adaptor, water bag, and coupler be disinfected before sterilization.

Country of		Chemical name/type	Ethanol
Country of origin :		Trade name	Ethanol
Any: Any c	,	Manufacturer	Any manufacturer
FRA: France	-	Country of origin	Any
DEU: Germ	,	Concentration	80%
	d Kingdom	(Dilution ratio)	(original solution)
USA: Unite	d States	Time	Wiping only
JPN: Japan		Temperature	Room temperature
AUS: Austr	alia	Humidity	Normal humidity
		Pressure	Normal pressure
Water bag/	Mode	el name <sup>*1</sup>	
coupler	UAWB-021A UAWB-022A		ОК
	UAGL-001AHA	UAGV-024A	
	UAGL-002AHA	UAGV-027A	
	UAGL004A	UAGV-028A	
	UAGL-004AHA	UAGV-029A UAGV-030A	—
Biopsy	UAGV020A	UAGV-030A UAGV-034A	
adaptor <sup>*2</sup> UAGV-023A UAGV-035A		UAGV-037A	
	UAGV-036A		—
		UAGV-031A	
	UAGL-015A	UAGV-032A	_
		UAGV-033A	
	UAFS-001A	UAFS-002A	
Magnetic	UAFS-003A	UAFS-004A	
sensor	UAFS-005A	UAFS-006A	_
securing	UAFS-007A UAFS-009A	UAFS-008A UAFS-010A	
adaptor	UAFS-009A UAFS-011A	UIFR-A500	

OK : Use of the chemical is permitted. \_\_\_\_ : Use of the chemical is not permitted.

Isopropyl alcohol				
Isopropyl alcohol				
Any manufacturer				
Any				
70%				
(original solution)				
Wiping only				
Room temperature				
Normal humidity				
Normal pressure				
ОК				
-				
_				
_				
-				

- \*1: Note that some of the listed models may not be available in some regions. For the transducers that support each item, refer to the operation manual supplied with the corresponding item.
- \*2: Sterilize the biopsy adaptor before and after use.

## **High-Level Disinfection -1 (Except for USA)**

\* Note that local regulations may require that the biopsy adaptor, water bag, and coupler be disinfected before sterilization.

Country of origin :		Chemical name/ type	Hydrogen peroxide, Peracetic acid	Chlorine	e dioxide
FRA: Franc	Any: Any country FRA: France DEU: Germany		PeraSafe <sup>™*3</sup>	Tristel Fuse for Instruments	Tristel Sporicidal Wipes/ Tristel Rinse Wipes
	ed Kingdom	Manufacturer	DuPont	Tristel Solut	ions Limited
USA: Unite		Country of origin	Any	G	BR
	JPN: Japan AUS: Australia		1.62% (61.7)	0.012% (1 sachet in 5 liters of water)	Working solution
		Time	10 min.	5 min.	Wiping only
		Temperature	Room temperature	Room ter	mperature
		Humidity	Normal humidity		humidity
		Pressure	Normal pressure	Normal	pressure
Water bag/ coupler UAWB-021A UAWB-022A		name <sup>*1</sup>	ОК	ОК	ОК
Biopsy adaptor <sup>*2</sup>	UAGL-001AHA UAGL-002AHA UAGL004A UAGL-004AHA UAGV020A UAGV-023A	UAGV-027A UAGV-028A	ОК	ОК	_
duptor	UAGV-035A UAGV-036A		ОК	ОК	_
	UAGL-015A	UAGV-031A UAGV-032A UAGV-033A	ОК	_	_
Magnetic sensor securing adaptor	UAFS-001A UAFS-003A UAFS-005A UAFS-007A UAFS-009A UAFS-011A	UAFS-002A UAFS-004A UAFS-006A UAFS-008A UAFS-010A UIFR-A500	ОК	ОК	_

OK : Use of the chemical is permitted. - : Use of the chemical is not permitted.

\*1: Note that some of the listed models may not be available in some regions. For the transducers that support each item, refer to the operation manual supplied with the corresponding item.

Except for	
Japan and China ,	
China / E	
日本・中国	
国以外向け /	
) / 日本・中国以外向け / 日本及中国以外的国家或地区适用	
<b>以外的国家</b> :	
或地区适用	

Succindialdehyde	Per	acetic acid	Peracetic acid, Hydrogen peroxide, Acetic acid
gigasept <sup>®</sup> FF (new)	NU-CIDEX <sup>®</sup>	Anioxyde 1000	gigasept <sup>®</sup> PAA concentrate
Schülke & Mayr GmbH	1&1	Laboratories ANIOS	Schülke & Mayr GmbH
DEU	GBR	FRA	DEU
5% (20)	0.35% (working solution)	0.15% (working solution)	2% (a double-chambered bottle into 5 liters of water)
15 min.	5 min.	30 min.	5 min.
Room temperature	Room	temperature	Room temperature
Normal humidity	Norn	nal humidity	Normal humidity
Normal pressure	Norr	nal pressure	Normal pressure
ОК	ОК	ОК	ОК
ОК	_	OK	ОК
ОК	_	ОК	ОК
ОК	_	OK	ОК
		_	OK

- \*2: Sterilize the biopsy adaptor before and after use.
- \*3: PeraSafe<sup>™</sup> may be classified as a sterilizing agent or as a high-level disinfectant depending on the local regulations.

# **High-Level Disinfection -2**

\* Note that local regulations may require that the biopsy adaptor, water bag, and coupler be disinfected before sterilization.

Country of	Country of origin :		Chemical name/type	Glutaraldehyde		
Any: Any country FRA: France		Trade name	CIDEX®	CIDEX PLUS <sup>®</sup> 28 day solution		
DEU: Germany GBR: United Kingdom USA: United States		Manufacturer	181	1&1		
			Country of origin		Any	
JPN: Japan AUS: Australia			Concentration (Dilution ratio)	2.4% (working solution)	3.4% (working solution)	
			Time	45 min.	20 min.	
			Temperature	Room temperature		
			Humidity	Normal humidity		
			Pressure	Normal pressure		
Water bag/		Model name <sup>*1</sup>				
coupler	UAWB-021A UAWB-022A			OK		
Biopsy adaptor <sup>*2</sup>	UAGL-001AHA UAGL-002AHA UAGL004A UAGL-004AHA UAGV020A UAGV-023A		UAGV-024A UAGV-027A UAGV-028A UAGV-029A		_	
			UAGV-030A UAGV-034A UAGV-037A			
	UAGV-035A UAGV-036A			_		
	UAGL-015A		UAGV-031A UAGV-032A UAGV-033A	_		
Magnetic sensor securing adaptor	sensor UAFS-005A securing UAFS-007A adaptor UAFS-009A		UAFS-002A UAFS-004A UAFS-006A UAFS-008A UAFS-010A		_	
	UAFS-011A		UIFR-A500			

OK : Use of the chemical is permitted. - : Use of the chemical is not permitted.

Glutaraldehyde					Ortho- phthalaldehyde
MetriCide™	MetriCide™ 28	MetriCide™ Plus 30	WAVICIDE <sup>®</sup> -01	Sporicidin <sup>®</sup>	CIDEX <sup>®</sup> OPA
Metrex Research, Inc			Medical Chemical Corporation	Contec Inc.	າສາ
	USA		USA	USA	Any
2.6% (working solution)	2.5% (working solution)	3.4% (working solution)	2.65% (working solution)	1.12% (working solution)	0.55% (working solution)
45 min.	90 min.	90 min.	45 min.	20 min.	12 min.
	25°C (77°F)		Room temperature	Room temperature	Room temperature
Ν	ormal humidi	ty	Normal humidity	Normal humidity	Normal humidity
N	ormal pressu	re	Normal pressure	Normal pressure	Normal pressure
ОК			ОК	ОК	ОК
_			_	_	ОК
_			_	_	ОК
_			_	_	ОК
_			_	_	OK

- \*1: Note that some of the listed models may not be available in some regions. For the transducers that support each item, refer to the operation manual supplied with the corresponding item.
- \*2: Sterilize the biopsy adaptor before and after use.

21

### Sterilization

\* Note that local regulations may require that the biopsy adaptor, water bag, and coupler be disinfected before sterilization.

Country of	Country of origin :		Chemical name/type	Ethylene oxide gas $^{*3}$
-	Country of origin :		Trade name	Ethylene oxide gas
	Any: Any country FRA: France		Manufacturer	Any manufacturer
			Country of origin	Any
DEU: Germ	-		Concentration	10%
	d Kingdom		(Dilution ratio)	1078
USA: Unite	d States		Time	Exposure time: 7 hours
JPN: Japan			_	Aeration time: 12 hours
AUS: Austr	alia		Temperature	50°C (122°F)
			Humidity	50%
			Pressure	980 hPa [gauge]
Water bag/		Model	name <sup>*1</sup>	
Water bag/ coupler	UAWB-021A			OK
coupier	UAWB-022A	١		
	UAGV-023A			ОК
	UAGL-001AHA		UAGV-024A	
	UAGL-002AHA		UAGV-027A	
	Biopsy JAGV020A		UAGV-028A	
			UAGV-029A	OK
Biopsy			UAGV-030A	
adaptor <sup>*2</sup>			UAGV-034A	
adaptor			UAGV-037A	
	UAGV-035A			_
	UAGV-036A			
			UAGV-031A	
	UAGL-015A		UAGV-032A	OK
			UAGV-033A	
	UAFS-001A		UAFS-002A	
Magnetic	UAFS-003A		UAFS-004A	
sensor	UAFS-005A		UAFS-006A	_
securing	UAFS-007A		UAFS-008A	
adaptor	UAFS-009A		UAFS-010A	
	UAFS-011A		UIFR-A500	

OK : Use of the chemical is permitted. - : Use of the chemical is not permitted.

- \*1: Note that some of the listed models may not be available in some regions. For the transducers that support each item, refer to the operation manual supplied with the corresponding item.
- \*2: Sterilize the biopsy adaptor before and after use.

Steam	Hydrogen peroxide plasma <sup>*4</sup>
High-pressure steam sterilization	STERRAD <sup>®</sup> 50/100S/200/NX
	۲&۱
Any	Any
	Use STERRAD system cassettes specifically designed for the sterilizer.
18 min.	Select [SHORT Cycle] for STERRAD 100S/200 <sup>*5</sup> and [STANDARD Cycle] for STERRAD NX.
134°C (273.2°F)	There are no sterilization related parameters
	(temperature, humidity, pressure, or time) to be
	set by the operator.
_	_
ОК	ОК <sup>*6</sup>
OK	ОК
ОК	ОК
ОК	ОК
_	ОК

- \*3: After sterilization, thoroughly degas the sterilized item to remove all gas residues.
- \*4: The screw on the biopsy adaptor may become tight after sterilization. This does not indicate that the biopsy adaptor has been damaged.
- \*5: For STERRAD 100S/200 for the USA, no Cycle option is provided.
- \*6: Use a "booster" for the lumen of the biopsy adaptor for STERRAD 50/100S/200.

# 4. Other Information

• For Europe

CANON MEDICAL SYSTEMS EUROPE B.V.

EC REP

Zilverstraat 1, 2718 RP Zoetermeer, The Netherlands

https://eu.medical.canon

For other countries, please contact your local distributor.

• Legal manufacturer



Canon Medical Systems Corporation 1385, Shimoishigami, Otawara-shi, Tochigi 324-8550, Japan